



## Test Mix - Engelhardt

### Conditions

Column: ACE 3 C18, 150 x 4.6mm  
Part Number: ACE-111-1546  
Mobile Phase: 55:45 MeOH/H<sub>2</sub>O  
Flow Rate: 1.0ml/min  
Temperature: Ambient  
Detection: UV, 254nm

### Compounds

1. Thiourea
2. Aniline
3. Phenol
4. Toluidine (o,m,p combined)
5. Dimethylaniline
6. Ethyl benzoate
7. Toluene
8. Ethylbenzene

