

Conditions

Column: ACE Excel 2 C18-PFP
Dimensions: 100 x 3.0 mm
Part Number: EXL-1010-1003U
Mobile Phase: A: 20 mM ammonium acetate pH 6.0 in H₂O
B: 20 mM ammonium acetate pH 6.0 in MeOH/H₂O (9:1 v/v)

Time (mins)	%B
0.0	6
5.5	63

Flow Rate: 1.2 mL/min
Injection: 2 µL
Temperature: 27 °C
Detection: UV, 220 nm

