

Conditions

Column: ACE 5 C18-PFP
Dimensions: 150 x 4.6 mm
Part Number: ACE-1210-1546
Mobile Phase: A: 20 mM Heptafluorobutyric acid in H₂O
B: 20 mM Heptafluorobutyric acid in MeCN

Application #AN6920

Time (mins)	%B
0	25
15	40
17	100
19	100
20	25
40	25

Flow Rate: 1 mL/min
Injection Volume: 30 µL
Temperature: 30 °C
ELSD: Evaporator Temp: 50 °C
Gain: 4
Filter: 10 s
Sampling rate: 30 Hz

1. Streptomycin
2. Amikacin
3. Tobramycin
4. Gentamicin (i, ii & iii)
5. Neomycin

