

## Conditions

Column:	ACE 10 C18
Dimensions:	250 x 4.6 mm
Part Number:	ACE-131-2546
Mobile Phase:	50 mM monobasic potassium phosphate pH 5.0: MeCN (24:1 v:v)
Flow Rate:	1.5 mL/min
Injection:	10 $\mu$ L
Temperature:	22 $^{\circ}$ C
Detection:	UV, 230 nm
Sample:	Amoxicillin Oral Suspension: 1 mg/mL amoxicillin in 50 mM monobasic potassium phosphate pH 5.0
System:	Chromaster 600

