

Analysis of Niuhuang Shangqing Pills

Application #AN6690

Conditions

Columns: ACE Excel 5 AQ
Dimensions: 250 x 4.6 mm
Part Numbers: EXL-126-2546U
Mobile Phase: A: 0.05 % phosphoric acid (aq.)
B: MeCN

Time (mins)	%B
0	10
18	23
30	27
35	35
40	35
45	50
50	10

Flow Rate: 1.0 mL/min
Injection: 10 µL
Temperature: 30 °C
Detection: UV, 240 nm
Method: ChP.2015, Vol. 1

