

Antihistamines and Expectorants – Mobile Phase Effects

ACE[®]
Ultra-inert
UHPLC & HPLC Columns

Application #AN3960

Conditions

Column: ACE 5 C18
Dimensions: 150 x 4.6 mm
Part Number: ACE-121-1546
Mobile Phase: (a) MeOH/50 mM KH₂PO₄ pH 3.0 in H₂O (45:55 v/v)
(b) MeCN/50 mM KH₂PO₄ pH 3.0 in H₂O (28:72 v/v)
Flow Rate: 1 mL/min
Temperature: Ambient
Detection: UV, 220 nm

