

BSA Tryptic Digest Profiling

Application #AN2000

Conditions

Column: ACE 5 C18-300
Dimensions: 150 x 4.6 mm
Part Number: ACE-221-1546
Mobile Phase: A: 1% TFA in H₂O
B: 1% TFA in MeCN/H₂O (1:1 v/v)

Time (mins)	%B
0	4
5	4
25	20
45	20
75	40
95	65
115	70
120	4

Flow Rate: 1 mL/min
Temperature: Ambient
Detection: UV, 214 nm

