

DNPH Derivatised Carbonyls (EPA Method TO-11A Test Mix)

Application #AN4700

Conditions

Column: ACE Excel 1.7 C18
Dimensions: 100 x 3.0 mm
Part Number: EXL-171-1003U
Mobile Phase: A: H₂O
B: MeCN

Time (mins)	%B
0.0	55
4.0	55
7.0	100
8.0	100
8.2	55
13.2	55

Flow Rate: 1 mL/min
Injection: 1 µL
Temperature: 50 °C
Detection: UV, 365 nm

