

Determination of Ferulic Acid in Ligusticum Chuanxiong

ACE[®]
Ultra-inert
UHPLC & HPLC Columns

Application #AN6120

Conditions

Column:	ACE UltraCore 2.5 SuperC18
Dimensions:	100 x 4.6 mm
Part Number:	CORE-25A-1046U
Mobile phase:	1% acetic acid (aq.):MeOH (70:30 v/v)
Flow Rate:	1.0 mL/min
Injection:	10 µL
Temperature:	Ambient
Detection:	UV, 321 nm
Monograph method:	ChP.2015, Vol. 1

