

Determination of Notoginseng Triterpenes in Xueshuantong Injection

Application #AN6490

Conditions

Column: ACE Excel 5 C18-AR
Dimensions: 150 x 4.6 mm
Part Number: EXL-129-1546U
Mobile Phase: A: H₂O
B: MeCN

Time (mins)	%B
0	20
45	20
65	34
85	34
86	90
96	90
97	20

Flow Rate: 1.0 mL/min
Injection: 10 µL
Temperature: 30 °C
Detection: UV, 203 nm
Method: National Medicine Standards WS3-B-3829-98-2011

1. Notoginsenoside R1
2. Ginsenoside Rg1
3. Ginsenoside Re
4. Ginsenoside Rb1
5. Ginsenoside Rd

