

Determination of Panax Notoginseng Saponins in Salvia Miltiorrhiza Tablets

Application #AN6450

Conditions

Column: ACE Excel 5 C18
Dimensions: 250 x 4.6 mm
Part Number: EXL-121-2546U
Mobile Phase: A: H₂O
B: MeCN

Time (mins)	%B
0	19
35	19
55	29
70	29
100	40

Flow Rate: 1 mL/min
Injection: 20 µL
Temperature: Ambient
Detection: UV, 203 nm
Method: ChP. 2015, Vol. 1

