

Conditions

Column: ACE Excel 1.7 C18-PFP
Dimensions: 2 x 100 x 3.0 mm (coupled)
Part Number: 2 x EXL-1710-1003U
Mobile Phase: A: 0.1% formic acid in H₂O
B: 0.1% formic acid in MeCN

Time (mins)	%B
0.00	5
0.72	5
15.72	50
18.72	100
20.72	100
22.72	5

Flow Rate: 0.8 mL/min
Injection: 2 µL
Temperature: 80 °C
Detection: UV, 254 nm

