

Herbicide Impurity Profile

Application #AN2130

Conditions

Column: ACE UltraCore 2.5 SuperC18
Dimensions: 150 x 4.6 mm
Part Number: CORE-25A-1546U
Mobile Phase: A: MeCN/H₂O/TFA (5:95:0.05 v/v/v)
B: MeCN/TFA (99.9:0.05 v/v)

Time (mins)	%B
0	10
3	10
35	100
55	100
56	10
60	10

Flow Rate: 0.6 mL/min
Injection: 10 µL
Temperature: 25 °C
Detection: UV, 240 nm

Technical Grade Herbicide

Active component

