

Lipid Classes Separation from *Drosophila Melanogaster*

Conditions

Column: ACE 3 SIL
Dimensions: 150 x 3.0 mm
Part Number: ACE-117-1503
Mobile Phase: A: IPA/MeCN (20:80 v/v)
B: IPA/0.02 M ammonium formate (20:80 v/v)

Time (mins)	%B
0.0	8
1.0	8
5.0	9
10.0	20
16.0	25
23.0	35
26.0	8

Flow Rate: 0.3 mL/min
Detection: LTQ Orbitrap MS
Positive and negative ion mode

