

NEW: **Ball-Valve Funnel with Hinged Lid**

Strong Protection when collecting
Liquid Waste!



Clever Technology, Strong Protection!

Safety Funnel



NEW: Improved Closing

Comfortable handling: easy, quick opening and closing with "click"-closure.



NEW: Hinged Lid with Sealing Ring

For secure closure, when the funnel is not in use.



NEW: Stable Hinge

A double-screw hinge, constructed of rust-free stainless steel: for long life and stability.



NEW: Optimized Spray Protection

Air ducts directed downwards ensure for secure ventilation and cause any spray created to be returned to the funnel.



NEW: Optimized Dirt Sieve

For collecting mixing magnets or large dirt particles. Laboratory bottles can also be placed here to dry off. The hole size has been optimized, in order to improve the running-off of chemicals. Both sieve and spray protection can be removed for cleaning.



Ball-Valve

The valve only opens up during filling, and then closes automatically afterwards - for a big "Plus" regarding safety!



Various Screw-Threads



Suitable for different sizes of collecting canister (see Table on Page 3).



Safety Lance



Ensures for safe and good flowing of sensitive chemicals, without any "stuttering."



NEW: Robust Design

The new design provides for a high degree of stability as well as long life, and prevents any deformation of the funnel from occurring.

Information for Ordering




Fig.	Part No.	Screw-Thread Of Container	Lance	Electrically conductive	Delivery with Adapter
-	MM317 642	GL 45	●	-	-
-	MM317 622	GL 45	●	●	-
-	MM317 638	GL 45	-	●	-
-	MM317 649	S 50	●	-	●
-	MM317 629	S 50	●	●	●
-	MM317 644	S 51	●	-	●
-	MM317 624	S 51	●	●	●
-	MM317 632	B 53	●	-	●
-	MM317 645	S 55	●	-	-
-	MM317 625	S 55	●	●	-
A	MM317 641	S 60/61	●	-	-
B	MM317 621	S 60/61	●	●	-
-	MM317 646	S 65	●	-	●
-	MM317 626	S 65	●	●	●
-	MM317 648	S 70/71	●	-	●
-	MM317 628	S 70/71	●	●	●
-	MM317 647	B 83	●	-	●
-	MM317 627	B 83	●	●	●
-	MM317 643	S 90	●	-	●
-	MM317 623	S 90	●	●	●
-	MM317 635	S 95	●	●	●
-	MM317 630	Dual Screw-Thread R2" BSP/ G2" + 2" Trisure	-	●	-
Models without ball valve, for highly viscous liquids					
-	MM317 633	GL 45	-	●	-
-	MM317 634	S 60/61	-	●	-



Information about "Delivery with Adapter:"

The safety funnel is available in various container screw-thread sizes. In order to be environmentally-friendly, certain variants are delivered together with a screw-thread adapter. Those to which this applies can be seen in the column "Delivery with Adapter."



* The total height can vary, dependent upon the thread adapter with which the funnel has been delivered.



Safety Funnel



**THE No. 1 System
for Solvent Safety**

MAC-MOD Analytical

103 Commons Court, Chadds Ford PA 19317



610-358-9696



efaden@mac-mod.com



mac-mod.com