

MAC-MOD Overview

Range of Canisters



Product Group 1: Electrically Conducting Canisters

Image	Thread	Liters	Material	Colour	UN-Certified	Special Feature	Inch W x H x D	Part Number
A	S 50	10	PE-HD-el	black	Yes	UN-X approval; with viewing strip	7.48" x 12.40" x 9.06"	MMSCCA2
B	S 60 / 61	10	PE-HD-el	black	Yes	UN-Y approval	7.28" x 10.43" x 11.41"	MMSCCA5
C	S 60 / 61	10	PE-HD-el	black	No	with float	7.28" x 10.43" x 11.41"	MMSCCA6
D	S 60 / 61	20	PE-HD-el	black	Yes	UN-Y approval	7.28" x 19.69" x 11.41"	MMSCCA12
E	S 60 / 61	20	PE-HD-el	black	No	with float	7.28" x 19.69" x 11.41"	MMSCCA13
F	S 60 / 61	30	PE-HD-el	black	Yes	UN-Y approval	9.44" x 17.91" x 14.33"	MMSCCA19
G	S 70 / 71	60	PE-HD-el	black	Yes	UN-Y approval	12.99" x 24.61" x 15.59"	MMSCCA27
H	S 70 / 71	60	PE-HD-el	black	No	with float	12.99" x 27.17" x 15.59"	MMSCCA26
I	S 90	10	PE-HD-el	black	Yes	UN-X approval; with viewing strip	7.68" x 14.96" x 7.68"	MMSCCA8

MAC-MOD Overview

Range of Canisters



Product Group 2: Non-Electrically Conducting Canisters

Image	Thread	Liters	Material	Colour	UN-Certified	Special Feature	Inch W x H x D	Part Number
A	GL 45	2.5	PE-HD	natural (white)	Yes	UN-Y approval	4.53" x 8.28" x 5.91"	MMSCCA10
B	GL 45	5	PE-HD	natural (white)	Yes	UN-Y approval	5.91" x 9.84" x 7.68"	MMSCCA20
C	GL 45	10	PE-HD	natural (white)	Yes	UN-Y approval	7.48" x 12.01" x 9.06"	MMSCCA1
D	S 51	5	PE-HD	natural (white)	Yes	UN-X approval	5.71" x 9.84" x 7.48"	MMSCCA21
E	S 55	5	PE-HD	natural (white)	Yes	UN-Y approval	6.30" x 9.06" x 7.28"	MMSCCA23
F	S 55	5	F-PE-HD	natural (white)	Yes	UN-X approval; with inner fluorination	5.71" x 9.84" x 7.48"	MMSCCA22
G	S 55	10	F-PE-HD	natural (white)	Yes	UN-X approval; with inner fluorination	7.48" x 12.20" x 9.06"	MMSCCA4
H	S 55	10	PE-HD	natural (white)	Yes	UN-Y approval	7.28" x 12.01" x 8.86"	MMSCCA3
I	S 60 / 61	12	PE-HD	natural (white)	Yes	UN-X approval	7.87" x 13.78" x 9.25"	MMSCCA9
J	S 60 / 61	20	PE-HD	natural (white)	Yes	UN-X approval	10.24" x 15.35" x 11.38"	MMSCCA11
K	S 60 / 61	20	PE-HD	natural (white)	No	with float	10.24" x 17.91" x 11.22"	MMSCCA14
L	S 60 / 61	30	PE-HD	natural (white)	Yes	UN-Y approval	11.42" x 15.75" x 14.96"	MMSCCA16
M	S 60 / 61 & GL 45	30	PE-HD	natural (white)	No	2 openings	11.42" x 15.75" x 14.96"	MMSCCA18
N	S 70 / 71	60	PE-HD	natural (white)	Yes	UN-Y approval	12.99" x 25.00" x 14.57"	MMSCCA25
O	S 90	10	PE-HD	natural (white)	Yes	UN-Y approval	7.68" x 14.96" x 7.68"	MMSCCA7
P	S 95	5	PE-HD	natural (white)	Yes	UN-Y approval	6.69" x 12.20" x 6.69"	MMSCCA24
Q	S 95	20	PE-HD	natural (white)	No		10.24" x 15.35" x 11.42"	MMSCCA15

MAC-MOD Overview

Range of Canisters



Product Group 5: Round Canisters

Image	Thread	Liters	Material	Colour	UN-Certified	Special Feature	Inch W x H x D	Part Number
A	B 53	2	PE-HD	natural (white)	No		4.69" x 10.24" x 4.69"	MMSCRC3
B	S 65	5	PE-HD	natural (white)	No	with handle and measuring scale	6.57" x 12.99" x 6.57"	MMSCRC6
C	S 65	10	PE-HD	natural (white)	No	with handle and measuring scale	8.07" x 16.93" x 8.07"	MMSCRC2
D	S 65	25	PE-HD	natural (white)	No	w. carrying handles and measuring scale	10.94" x 22.83" x 10.94"	MMSCRC4
-	S 65	60	PE-HD	natural (white)	No	w. carrying handles and measuring scale	13.78" x 32.48" x 13.78"	MMSCRC8
E	B 83	4	PE-HD	natural (white)	No		6.10" x 13.31" x 6.10"	MMSCRC5
F	B 83	10	PP	natural (white)	No	with carrying handles	9.84" x 15.35" x 9.84"	MMSCRC1
G	B 83	50	PP	natural (white)	No	with carrying handles	14.96" x 26.77" x 14.96"	MMSCRC7

Product Group 6: Flat Canisters

Image	Thread	Liters	Material	Colour	UN-Certified	Special Feature	Inch W x H x D	Part Number
H	B 63	8	PP	natural (white)	No	with drain valve	14.76" x 6.50" x 12.01"	MMSCFC2
I	S 70 / 71	20	PE-HD	natural (white)	Yes	UN-Y approval	11.14" x 11.50" x 14.80"	MMSCFC1