

Non-Steroidal Anti-Inflammatory Drugs – Fast Analysis

ACE[®]
Ultra-inert
UHPLC & HPLC Columns

Application #AN2080

Conditions

Column: ACE Excel 2 SuperC18
Dimensions: 50 x 3.0 mm
Part Number: EXL-1011-0503U
Mobile Phase: A: 0.1% formic acid in H₂O
B: 0.1% formic acid in MeCN

| Time (mins) | %B |
|-------------|----|
| 0.00 | 20 |
| 0.25 | 20 |
| 3.50 | 70 |
| 4.00 | 20 |
| 4.25 | 20 |

Flow Rate: 0.86 mL/min
Injection: 1.4 µL
Temperature: 40 °C
Detection: UV, 254 nm

