

Pharmaceutically Relevant Compounds (V)

ACE[®]
Ultra-inert
UHPLC & HPLC Columns

Application #AN2500

Conditions

Column: ACE 3 C18-PFP
Dimensions: 50 x 2.1 mm
Part Number: ACE-1110-0502
Mobile Phase: A: 0.1% formic acid in H₂O
B: 0.1% formic acid in MeOH

Time (mins)	%B
0	3
5	100

Flow Rate: 0.6 mL/min
Temperature: 40 °C
Detection: UV, 254 nm

