

Vitamin B Mixture in HILIC

Application #AN5050

Conditions

Column: ACE 1.7 HILIC-N
Dimensions: 100 x 2.1 mm
Part Number: HILN-17-1002U
Mobile Phase: A: 10 mM ammonium formate pH 3.0 in MeCN/H₂O (50:50)
B: 10 mM ammonium formate pH 3.0 in MeCN/H₂O (90:10)

Time (mins)	%B
0	100
4	100
15	0
16	0
17	100
37	100

Flow Rate: 0.21 mL/min
Injection: 2 µL
Temperature: 22 °C
Detection: UV, 260 nm

