COMPREHENSIVE CONVERSION TABLES TO EXISTING CHIRAL MANUFACTURERS

MAC-MOD Chiral	Daicel/Chiral Technologies	Phenomenex	Regis	Dr. Maisch	ES Industries
Amy-1	CHIRALPAK® AD	Lux Amylose-1	REFLECT C-Amylose A	ReproSil Chiral-AM	ChromegaChiral CCA
Amy-2	CHIRALPAK® AY	Lux Amylose-2	-	ReproSil Chiral-YM	ChromegaChiral CC3
Amy-3	CHIRALPAK® AS	-	-	ReproSil Chiral-AMS	ChromegaChiral CCS
Amy-4	-	-	-	-	-
Amy-5	CHIRALPAK® AZ	-	-	ReproSil Chiral-ZA	-
Amy-6	-	-	-	-	-
Amy-8	-	-	-	-	-
Cel-1	CHIRALCEL® OD	Lux-Cellulose-1	REFLECT C-Cellulose B	ReproSil Chiral-OM	ChromegaChiral CCO
Cel-2	CHIRALCEL® OZ	Lux-Cellulose-2	-	ReproSil Chiral-ZM	ChromegaChiral CC2
Cel-3	CHIRALCEL® OJ	Lux-Cellulose-3	-	ReproSil Chiral-JM	ChromegaChiral CCJ
Cel-4	CHIRALCEL® OX	Lux-Cellulose-4	-	ReproSil Chiral-XM	ChromegaChiral CC4
Cel-9	-	-	-	-	-
Amy-4F	-	-	-	-	ChromegaChiral CCA F4
Cel-2F	-	-	-	-	ChromegaChiral CCO F2
Cel-4F	-	-	-	-	ChromegaChiral CCO F4
Cel-FT	-	-	-	-	ChromegaChiral CCO F4 T3
Amy -1-3B	-	-	-	-	ChromegaChiral CCX
Amy – 3-5B	-	-	-	-	ChromegaChiral CCU
Cel – 2,Cl-E	-	-	-	-	ChromegaChiral CCC
Im-Amy-1	CHIRALPAK® IA	Lux-i-Amylose-1	REFLECT I-Amylose-A	ReproSil Chiral -MIA	-
Im-Amy-3	CHIRALPAK® IH	-	-	-	-
Im-Amy-5	CHIRALPAK® IF			ReproSil Chiral -MIF	
Im-Amy-6	CHIRALPAK® IG	Lux-i-Amylose-3	-	ReproSil Chiral-MIG	-
Im-Cel-1	CHIRALPAK® IB	-	REFLECT I-Cellulose-B	ReproSil Chiral -MIB	-
Im-Cel-C	CHIRALPAK® IC	Lux-i-Cellulose-5	REFLECT I-Cellulose-C	ReproSil Chiral -MIC	-
Im-Cel-2	CHIRALPAK® IM	-	REFLECT I-Cellulose-Z		-
Im-Cel-9	CHIRALPAK® IK	-	-	-	-

CHIRALCEL® & CHIRALPAK® are registered trademarks of DAICEL Chemical Industries, Ltd. of Japan. AD, AY, AS, OD, OZ, OJ, OX, IA, IH, IF, IB, IC, IK are registered trademarks of DAICEL Chemical Industries, Ltd. of Japan. Lux is a registered trademark of Phenomenex, INC. REFLECT is a registered trademark of Regis Technologies. ReproSil is a registered trademark of Dr. Maisch and Chromegachiral is a registered trademark of PERKIN-ELMER. Disclaimer: Comparative and representative separations may not be representative of all applications.